



meat selection

BEWI DOG® Meat Selection contains juicy meat with a firm texture. Just the way dogs like it. The wide range of flavours provides a welcome variety.



rich in tripe

Composition: 70 % Meat and innards: Tripe (15 %), Meat, Lung, Liver, Udder; Salmon oil (0.5 %); Egg shells, dried; Minerals

Jar 400g Pack of six Article No. 518415
Jar 800g Pack of six Article No. 518815

Analytical constituents:

Protein	10.5 %	Fat content	5.5 %
Crude fibre	0.4 %	Crude ash	2.2 %
Moisture	79 %		



rich in lamb

Composition: 70 % Meat and innards: Meat (15 %), Lamb offals, Heart, Liver, Kidney, Udder; Salmon oil (0.5 %); Egg shells, dried; Minerals

Jar 400g Pack of six Article No. 518455
Jar 800g Pack of six Article No. 518855

Analytical constituents:

Protein	10.5 %	Fat content	5.5 %
Crude fibre	0.4 %	Crude ash	2.0 %
Moisture	79 %		



rich in chicken

Composition: 70 % Meat and innards: Chicken meat (15 %), Lung, Liver, Kidney, Udder, Heart; Salmon oil (0.5 %); Egg shells, dried; Minerals

Jar 400g Pack of six Article No. 518425
Jar 800g Pack of six Article No. 518825

Analytical constituents:

Protein	10.5 %	Fat content	5.5 %
Crude fibre	0.4 %	Crude ash	2.2 %
Moisture	79 %		



rich in chicken

Composition: 70 % Meat and innards: meat, Beef lung (15 %), Liver, Kidney, Udder, Heart; Salmon oil (0.5 %); Egg shells, dried; Minerals

Jar 400g Pack of six Article No. 518465
Jar 800g Pack of six Article No. 518865

Analytical constituents:

Protein	11.0 %	Fat content	6.5 %
Crude fibre	0.3 %	Crude ash	2.0 %
Moisture	79 %		



meat selection



rich in venison

Composition: 70 % Meat and innards: Venison meat (15 %), Lung, Liver, Heart, Kidney, Udder; Salmon oil (0.5 %); Egg shells, dried; Minerals

Jar 400g Pack of six Article No. 518435
Jar 800g Pack of six Article No. 518835

Analytical constituents:

Protein	10.5 %	Fat content	5.5 %
Crude fibre	0.4 %	Crude ash	2.1 %
Moisture	79 %		



rich in veal

Composition: 70 % Meat and innards: Meat (15 %), Calf's lung, Heart, Liver, Kidney, Udder; Salmon oil (0.5 %); Egg shells, dried; Minerals

Jar 400g Pack of six Article No. 518445
Jar 800g Pack of six Article No. 518845

Analytical constituents:

Protein	10.5 %	Fat content	5.5 %
Crude fibre	0.4 %	Crude ash	2.2 %
Moisture	79 %		



rich in poultry

Composition: 70 % Meat and innards: Poultry meat (15 %), Liver, Lung, Heart, Kidney, Udder; Salmon oil (0.5 %); Egg shells, dried; Minerals

Jar 400g Pack of six Article No. 518405
Jar 800g Pack of six Article No. 518805

Analytical constituents:

Protein	10.5 %	Fat content	5.5 %
Crude fibre	0.4 %	Crude ash	2.0 %
Moisture	79 %		



rich in poultry hearts

Composition: 70 % Meat and innards: Poultry heart, uncut (15 %), Lung, Liver, Kidney, Rind, Udder; Salmon oil (0.5 %); Egg shells, dried; Minerals

Jar 400g Pack of six Article No. 518475
Jar 800g Pack of six Article No. 518875

Analytical constituents:

Protein	11.0 %	Fat content	6.5 %
Crude fibre	0.3 %	Crude ash	2.0 %
Moisture	79 %		